



# FOR YOUR CHILDREN

## Preventive dental care in 3 easy steps

We can make an accurate, quick and painless diagnosis on the dental health of your children, discovering the areas at highest risk of tooth decay



### DIAGNOcam - "Caries detection without X-rays at your daily practice"

Do you know is there a new radiation free diagnostic device for the caries detection that gives the same results of a device that uses Xrays? DIAGNOcam allows to have more certain diagnosis without the use of Xrays.

DIAGNOACAM uses the structures of the tooth as light conductors. Simultaneously, a digital video camera, captures the context.



### Oral hygiene control

After analyzing the teeth using Diagnokam device, we can move to the control of the oral hygiene techniques.

Passed also this second step now we move to the next one...



### Tri Plaque ID Ge

Three tone plaque disclosing gel can identify new mature and acid producing biofilms.

GC Tri Plaque ID Gel can not only differentiate between old and new plaque in a few easy steps, this unique gel can also highlight exactly where the bacteria are most active by disclosing the acidic pH.

Plaque is revealed on teeth in three tones:

- red/pink - newly formed plaque
- blue/purple - old plaque (more than 48 hours)
- light blue - high risk plaque



*Studio dentistico De Pascalis*

Via Burlamacchi, n. 4 - Milano - [www.studiodepascalis.eu](http://www.studiodepascalis.eu)